

VINYL CORPAK™ COLD THERAPY PACKS

- Intense cold therapy eases pain and helps reduce inflammation
- Remains flexible, comfortably contours to painful areas
- Durable, latex free exterior easily wipes clean
- Nontoxic, biodegradable fill maintains therapeutic temperature for 20 minutes
- Various sizes can accommodate any treatment area
- Store in freezer for convenience
- Customized packs can be used as a promotional item (call for details)



Available From:



INTENSE COLD THERAPY

The Vinyl CorPak™ effectively provides intense cold therapy to help relieve headaches, ease pain, and help reduce inflammation associated with strains, sprains, bruises and other soft tissue injuries. Conveniently store packs in the freezer for cold therapy the moment you need it! These packs stay soft and flexible, even when frozen, to comfortably contour to areas of discomfort for targeted cold therapy.

PROVIDES SOOTHING RELIEF

Help reduce inflammation and swelling by using the Vinyl CorPak to provide soothing relief from muscle and joint pain, back, shoulder, neck, and knee discomfort. Use for sports injuries, post-workout recovery, and rehabilitation. This pack is also ideal for fever reduction, chronic pain, bursitis, tendinitis, and arthritis.

IDEAL FOR CLINICAL AND HOME USE

Our packs contain more fill than other conventional ice packs, giving them the powerful ability to retain therapeutic temperatures longer. The nontoxic, biodegradable fill remains pliable when frozen, unlike most other cold therapy packs, while helping reduce environmental impact.

The durable, professional grade, latex free exterior easily wipes clean, making it ideal for clinical and home use. The smooth vinyl cover resists tearing and splitting, while the durable RF welded seam provides a secure seal to help prevent leaks.

SPECIFICATIONS

Product Selection: Vinyl CorPak™ Cold Therapy Packs

Materials: Food grade propylene glycol USP, Bentonite clay, Water
Cover: polyether urethane film

Latex: This product is not made with natural rubber latex

Country of Origin: Made in USA

California Residents Only:
⚠️ **WARNING:** Cancer and Reproductive Harm – www.P65Warnings.ca.gov

VINYL CORPAK™ COLD THERAPY PACKS

INSTRUCTIONS FOR USE

Simply store the pack in the freezer for a minimum of 1 hour prior to use.

Wrap pack in fabric cover before applying to treatment area.

The pack will maintain therapeutic temperature for 20 minutes.

All Vinyl CorPak's are:



New!

Pictured Above:

X-Large Vinyl CorPak™ (12.5" x 18.5")

Available From:



Medium

ACC-565
Product Dimensions: 6" x 10"
Product Weight: 1.5 lbs
Case Quantity: 16
Case Weight: 27 lbs
Case Dimensions: 12" x 12" x 6"



Large

ACC-566
Product Dimensions: 10" x 13"
Product Weight: 3.5 lbs
Case Quantity: 6
Case Weight: 30 lbs
Case Dimensions: 12" x 12" x 6"



Cervical

ACC-567
Product Dimensions: 6" x 20"
Product Weight: 2.5 lbs
Case Quantity: 10
Case Weight: 32 lbs
Case Dimensions: 12" x 12" x 6"



Tri-Sectional

ACC-568
Product Dimensions: 11" x 15"
Product Weight: 3.75 lbs
Case Quantity: 6
Case Weight: 25 lbs
Case Dimensions: 12" x 12" x 6"



X-Large

ACC-569
Product Dimensions: 12.5" x 18.5"
Product Weight: 5.8 lbs
Case Quantity: 4
Case Weight: 25 lbs
Case Dimensions: 12" x 12" x 6"