

# Easy as **1 2 3** CaviWipes™

Fast and Effective Cleaning for Safer Client Care



## CaviWipes HP

### Features (S.A.F.E.)

#### Simple

1-minute universal contact time.

#### Alcohol-free

Hydrogen peroxide suitable for sensitive patient populations.

#### Fast

1-step cleaning/disinfection.<sup>1</sup>

#### Effective

Disinfection against 67 pathogens including multidrug-resistant organisms (MDRO).

#### Part No. Description

16-1100	CaviWipes HP (6" x 6.75") 160 towelettes per canister
16-1150	CaviWipes HP XL (9" x 12") 65 towelettes per canister

## CaviWipes 2.0

### Features

#### Simple

2-minute universal contact time.

#### Agreeable

Neutral pH, low alcohol/quat formulation agreeable with healthcare surfaces and materials.

#### Fast

1-step cleaning/disinfection.<sup>1</sup>

#### Effective

Disinfection against 60 pathogens including MDRO.

#### Part No. Description

14-1100	CaviWipes 2.0 (6" x 6.75") 160 towelettes per canister
14-1150	CaviWipes 2.0 XL (9" x 12") 65 towelettes per canister

## CaviWipes Bleach

### Features

#### Simple

3-minute universal contact time.

#### Alcohol-free

1:10 bleach dilution ratio recommended by the CDC.

#### Fast

1-step cleaning/disinfection.<sup>1</sup>

#### Effective

Disinfection against 44 pathogens including *C. diff* spores and Canine Parvovirus.

#### Part No. Description

13-9100	CaviWipes Bleach (6" x 10.5") 90 towelettes per canister
13-1175	CaviWipes Bracket Fits all CaviWipes canisters

To learn more scan the QR code or visit [metrex.com](http://metrex.com)



1. Always follow product label for special precleaning and PPE instructions.